

### Cisa 1 Electromagnetic and digital sieve shaker BA 200N



**Note:**

Other units for different sieve diameters or amplitude of vibration available on request

Electromagnetic and digital sieve shaker for test sieves of 200mm.Ø/ 203.Ø.

Precision and reliability

The latest technology with intermittent sieving.

Includes closing quick system and safe system "cisa Press".

Recommended for materials with density >0,8 gr/cm

Suitable for both dry and wet sieving kit.

Movement induced electromagnetically, combining a Rotatory displacement with a vertical one of the Sample.

Maximum amplitude of vibration 2,00 mm.

Dimensions:400x310x95 (without sieves).

Net weight: 27 kg

Order No.	Description
00TZBA200N	Electromagnetic and digital sieve shaker BA 200N

### Cisa 2 Electromagnetic and digital sieve shaker BA 300N



Electromagnetic sieve shaker for test sieves of 300mm.Ø/200mm.Ø

Designed to separate solid and thick materials.

Includes the closing quick system "cisa press".

Sieve shaker with 3D movement.

The latest technology with intermittent sieving.

Precision and reliability.

No maintenance.

Strong, reliable and polyvalent.

Maximum amplitude of vibration 2,00 mm.

Four programs, for different kinds of sieving.

Noise level (with 4 sieves of 300mmØ and 500 gr. of sample)

Dimensions: 445x 335x 150mm.

Net weight 95 kg.

Order No.	Description
00TZBA300N	Electromagnetic and digital sieve shaker BA 300N

# Chromatography Consulting Service:

[consultas@scharlab.com](mailto:consultas@scharlab.com)

