



Ultrasonic cleaning baths "Ultrasons" without heating



CONTROL PANEL

Timer switch for 0 -15 minutes or continuous operation. Indicator lamp when the unit is in operation.

MODELS

Part No.	Capacity litres	Height / Width / Depth (usable) cm	Height / Width / Depth (exterior) cm	Generator power W	Timer 0-15'	Waste drain	Weight Kg
3000512 1	1	6 14 14	20 20 18	50	Yes	No	2,5
3000683 2	2,6	9 23 13	22 29 19	100	Yes	No	4
3001208 8	4	10 30 14	24 33 18	100	Yes	Yes	7
3000513 3	6	15 30 14	30 36 20	150	Yes	Yes	9,5
3000837 5	6,5	14 24 Ø	30 27,5 Ø	150	Yes	No	8
3000617 4	9	15 50 14	30 56 20	200	Yes	Yes	10,7
3000514 6	9	15 30 23	30 36 30	200	Yes	Yes	13
3000515 7	20	15 50 30	30 56 36	400	Yes	Yes	22

COMPLEMENT



Detergent "BIO-SEL"

Bio-degradable, phosphate free. Excellent cleaning power just by mixing only 20ml per litre of water. BIO-SEL, in a 2% solution has a 9.7 pH level. Low toxicity. This detergent disperses, emulsifies, segregates, dissolves, suspends and breaks down dirtiness.

- 1000025 Capacity 1 litre.
- 1000026 Capacity 4 litres.
- 1000852 Box of 12 x 1 litre bottles.
- 1000853 Box of 4 x 4 litres bottles.

Baskets

Lids are not compatible with baskets.

Lid with concentric reduction ring

ACCESSORIES

Baskets and Lids.
Made of AISI 310 stainless steel.

Cleaner model Part No.	Part No	BASKETS Height / Width / Depth (cm)	LID WITH CONCENTRIC REDUCTION RING Part No.	Number of positions	LID WITH HANDLE Part No.
3000512	6005123	5 11 11	6005122	1	6005121
3000683	6006833	7 20 12	6006832	1	6006831
3001208	6001209	7 26 13	-	-	6005131
3000513	6005133	12 26 13	6005132	2	6005131
3000837	6008373	13 23 Ø	-	-	6008371
3000617	6005144	12 46 12	6005134	3	6005124
3000514	6005143	13 29 23	6005142	2	6005141
3000515	6005153	13 49 29	6005152	6	6005151



Ultrasonic Bath "Ultrasons-Digit"

5 CYCLES ELECTRONIC DIGITAL TIME CONTROL
CAPACITY 1.4 LITRES

FEATURES

Internal tank made from stainless steel, with minimum and maximum level indicator.
External case, basket and lid made from PVC.
Maximum cleaning power, frequency 42 KHz

CONTROL PANEL

Pushbutton start/ stop
Timer 5 programs selection in seconds:
90" - 180" - 280" - 380" and 480".
Wide application range from jewellery, watches, optical components. Dentistry and general laboratory use.

MODEL

Part No.	Capacity litres	Height / Width / Depth (usable) cm	Height / Width / Depth (exterior) cm	Power W	Weight Kg
5900650	1.4	6 13 16	15 17 22	70	1.75

Comes complete with basket and see through transparent PVC lid.



Model "Elegant".