

043-100176 Vortex kémcsőrázó (vált. sebességű) ZX3

The ZX3 is a basic vortex mixer with adjustable stirring speed and two operating modes.

Electronic speed regulation: from 50 to 2400 rpm

Operating modes: touch, continuous

GENERAL FEATURES

Construction material: epoxy painted aluminum structure

Support system: 4 holdfast feet for strong fixing and high stability

Orbitale diameter: 5 mm

Power: 45 W

Power supply: 230 V / 50-60 Hz

Weight: 2 Kg

Dimension (WxHxD): 150x134x150 mm

