



### Polypropylene white coat with Mao collar

Standard colour: white, made of polypropylene, 35g/m<sup>2</sup> density, elasticated cuffs, velcro fastening, and Mao collar.

Description	Size	Pack (u.)	Art. No.
No pockets	XL	50	1115-22144

## CLOTHING: Workwear



### Certified laboratory white coat

Certified protection white coat.  
 It is made of 35% cotton, 65% polyester, with a PTFE coating, and has a lapel.  
 Unisex model.  
 Impermeable to particles, spraying and splashes of limited intensity.  
 Complies with norms UNE EN 340.2004, UNE EN 13034.2005 for PB chemical clothes.

#### With open cuff

Size	Pack (u.)	Art. No.
M	1	00000019-M
L	1	00000019-L
XL	1	00000019-XL



#### With elastic cuff

Size	Pack (u.)	Art. No.
M	1	00000020-M
L	1	00000020-L
XL	1	00000020-XL



### Laboratory white coat

The laboratory coat provides a protection barrier against liquids, biological substances and splashes of weak chemical products.  
 It has long-sleeves and is made of 65% polyester and 35% cotton.  
 With two lower pockets and one upper pocket, and a lapel.  
 It has a belt on the back and closes with buttons.

#### Lab coat for men

Size	Pack (u.)	Art. No.
46	1	1010070046
48	1	1010070048
50	1	1010070050
52	1	1010070052

#### Lab coat for ladies

Size	Pack (u.)	Art. No.
42	1	1010090842
44	1	1010090844
46	1	1010090846
48	1	1010090848