

[Click Below](#)

For General Use

Manual pipettes
Pip. Filler/ Stepper
Electronic Dispenser
Magnetic Stirrers
Shakers

For Analytical Labs

Electronic Titrator
Manual Dispenser
Rotary Evaporator
Hotplate
Overhead Stirrers

For Life Science Labs

Electronic Pipettes
SafeVac
Centrifuges-LS
Mixers / Vortex
Thermo Mixers

For Clinical Labs

Centrifuges-Clinical
Rollers and Rotators
Thermo Cyclers
Spectrophotometers



Info@dlabsci.com

High Speed Micro Centrifuges

D1008/D1008E/D2012

DLAB
Your Accessible Solution

**Little Space
and Great
Price**



D1008

Quick centrifugation in small space

For frequently useful tube quantities

For frequently used velocities

Easy to use

Quiet and consumes very low energy

Colors your lab

Small footprint

For strips and 8x2ml tubes (2rotors)

5000/ 7000 rpm (2models)

Start/stop

< 46 dB

Available in 4 colors with a great price advantage

D2012 Plus

Powerful centrifugal force 15100xg

Max speed 15000rpm

For upto 12 x 2/1.5/0.5/0.2ml tubes at 15,000 rpm

Single rotor for different capacity tubes

Refrigerated version available



[Click Below](#)

For General Use

Manual pipettes
Pip. Filler/ Stepper
Electronic Dispenser
Magnetic Stirrers
Shakers

For Analytical Labs

Electronic Titrator
Manual Dispenser
Rotary Evaporator
Hotplate
Overhead Stirrers

For Life Science Labs

Electronic Pipettes
SafeVac
Centrifuges-LS
Mixers / Vortex
Thermo Mixers

For Clinical Labs

Centrifuges-Clinical
Rollers and Rotators
Thermo Cyclers
Spectrophotometers



Web Page



Info@dlabsci.com

High Speed Micro Centrifuges

D1008/D1008E/D2012

DLAB
Your Accessible Solution

D1008

Technical data

Max. Speed	7000rpm
Max. RCF	2680×g
Rotor Capacity	0.2/0.5/1.5/2.0mL×8; 0.2mL×16 PCR strips or 0.2mL×2 PCR 8 strips
Run Time	Continuous operation
Driving Motor	DC motor
Power	AC110-240V/50Hz/60Hz
Noise level	≈45dB
Dimensions [D × W × H]	160 × 170 × 122mm
Weight	0.5kg

D1008E

Technical data

Max. Speed	5000rpm
Max. RCF	1360×g
Rotor Capacity	0.2/0.5/1.5/2.0mL×8; 0.2mL×16 PCR strips or 0.2mL×2 PCR 8 strips
Run Time	Continuous operation
Driving Motor	DC motor
Power	AC110-240V/50Hz/60Hz
Noise level	≈45dB
Dimensions [D × W × H]	160 × 170 × 122mm
Weight	0.5kg
Certification	CE cTÜVus FCC MCA

D2012

Technical data

Max. Speed	15000rpm (500-15000rpm), increment: 100rpm
Max. RCF	15100×g, increment: 100×g
Speed Accuracy	±20 rpm
Rotor Capacity	0.2mL/0.5mL/1.5mL/ 2mL×12
Run Time	30sec~99min/Continuous
Driving Motor	Brushless DC motor
Safety Devices	Door interlock, Over-speed detection, Over-temperature detection, Automatic internal diagnosis
Power	Single-phase, AC110V-240V, 50Hz/60Hz, 3A
Noise level	≈56dB
Acceleration/Braking time	25s?/ 25s?
Dimensions [D × W × H]	255 × 245 × 140mm
Weight	6kg
Additional Features	Speed/RCF switch, Short-time run function, Sound-alert function

Click Below

For General Use

Manual pipettes
Pip. Filler/ Stepper
Electronic Dispenser
Magnetic Stirrers
Shakers

For Analytical Labs

Electronic Titrator
Manual Dispenser
Rotary Evaporator
Hotplate
Overhead Stirrers

For Life Science Labs

Electronic Pipettes
SafeVac
Centrifuges-LS
Mixers / Vortex
Thermo Mixers

For Clinical Labs

Centrifuges-Clinical
Rollers and Rotators
Thermo Cyclers
Spectrophotometers



Web Page



Info@dlabsci.com

High Speed Micro Centrifuges

D1008/D1008E/D2012



Your Accessible Solution

Ordering Information

☆ Centrifuge - D1008 *EZeeMini* Centrifuge Always mention the needed rpms

Cat. No.	Descriptions
9031001012	Mini Centrifuge with DragonLab blue lid, including two rotors of A08-2 & A02-PCR8
9031002012	Mini Centrifuge with green lid, including two rotors of A08-2 & A02-PCR8
9031003012	Mini Centrifuge with pink lid, including two rotors of A08-2 & A02-PCR8
9031004012	Mini Centrifuge with yellow lid, including two rotors of A08-2 & A02-PCR8

Accessories of D1008

Cat. No.	Descriptions
19400015	Plastic rotor (blank LOGO), 7000rpm, 2ml/1.5ml*8, used with
19400017	Plastic rotor (blank LOGO), 7000rpm, PCR8*2, used with D1008
19400029	2ml to 0.2ml rotor adapter, used with D1008 rotor, 8pcs/pk
19400026	2ml to 0.5ml rotor adapter, used with D1008 rotor, 8pcs/pk

☆ Centrifuge - D2012 High Speed Mini Centrifuge

Cat. No.	Descriptions
9032002121	High Speed Mini Centrifuge (CE marked), with A12-2P plastic rotor kit

Accessories of D2012

Cat. No.	Descriptions
19400010	Plastic rotor kit with lid (blank LOGO), fixing clips and O'ring, 15000rpm, 2ml*12, used with D2012 plus
19500001	2ml to 0.2ml rotor adapter, used with A12-2, A12-2P, AS24-2 and A24-2P rotors, 24pcs/pk
19500002	2ml to 0.5ml rotor adapter, used with A12-2, A12-2P, AS24-2 and A24-2P rotors, 24pcs/pk