



## LABORATORY SCISSORS

Made of AISI 410 stainless steel.

### 1 Scissors universal "Multi purpose" with a serrated/notched blade.

Suitable for cutting hard plastic, cardboard and soft wire etc.

Part No. **1001468** Cut length 50 mm. Total length 180 mm.

### 2 Large laboratory scissors.

Part No. **1001469** Cut length 100 mm. Total length 250 mm.

### 3 Multipurpose straight scissors, general use and bandages.

Part No. **1001470** Cut length 50 mm. Total length 140 mm.

### 4 Multipurpose scissors, general use and bandages.

Part No. **1000806** Cut length 40 mm. Total length 145 mm.

### 5 Dissection scissors.

Part No. **1000807** Cut length 35 mm. Total length 125 mm.

### 6 Mini dissection scissors.

Part No. **1001471** Cut length 25 mm. Total length 110 mm.

### 7 Flat scissors with serrated edge, ideal for laminates, negatives and chromatography.

Part No. **1001472** Flat length 40 mm. Total length 160 mm.



## DISSECTION NEEDLES

Plastic handled, AISI 304 stainless steel needle.

Part No.	Point	Total length
<b>1000431</b>	Lancet	140 mm
<b>1000275</b>	Fine	140 mm



## DISSECTION NEEDLES

Hexagonal metallic handle with a AISI 304 stainless steel needle.

Part No.	Point	Total length
<b>1001359</b>	Lancet	40 mm
<b>1001360</b>	Fine	140 mm



## SCALPEL HANDLE No. 4

Made of AISI 304 stainless steel, comes with three sterile blades No. 20, 22 and 23.

Part No. **1001361**  
 Spares:  
 Part No. **1001362** 10 sterile blades n° 20  
**1001363** 10 sterile blades n° 22  
**1001364** 10 sterile blades n° 23



## ANATOMY SCALPEL

Made of AISI 304 stainless steel.

Blade length: 35 mm.

Part No. **1000809**

## FORCEPS RING SHAPED TIP

AISI 304 stainless steel forceps for picking up swabs and textiles etc.

Part No. **1001192**



## DISSECTION KIT

Model Student.

Contains: 1 fine point probe, 1 lancet probe, 1 blunt end forceps, 1 fine tip forceps, 1 scalpel handle n° 4, 3 sterile blades for scalpel handle n° 4, 1 scissors, 1 magnifying glass, 1 style probe, 1 scalpel.

Part No. **1001365**



## DISSECTION KIT

Wooden box.

Contents: 1 fine needle, 1 lancet needle, 1 dissection forceps blunt end, 1 forceps fine point, 1 scalpel handle n° 4 with 3 sterile blades Nos. 20, 22 & 23, 1 scalpel and 1 scissors.

Part No. **1001152**