



### DISSECTING FORCEPS AND GENERAL LABORATORY USE

Made of AISI 304 stainless steel with electro-bright finish.

Type	Description	Total length	Part No.
1	Dissecting Straight forceps round tips	140 mm	1000580
2	Dissecting Straight forceps fine tips	140 mm	1000801
3	Dissecting curved forceps round tips	140 mm	1000865
4	Dissecting curved forceps fine tips	140 mm	1000864

5	Sprung straight forceps with curved fine tips with wood laminate sides	160 mm	1001190
6	Forceps straight serrated fine point, with wood laminate sides.	160 mm	1001189
7	Laboratory forceps curved point	160 mm	6000128
8	Laboratory forceps fine point	120 mm	1001191



### KOLLE HANDLES


Nickel metal handle with loop end of Ni-Cr.

Handle length: 43 mm. Internal loop Ø: 3 mm.

Part No.

**1000033** Bakelite handle. Handle length: 200 mm.

**1000488** "Pasteur", with metal handle. Handle length: 290 mm.



#### 1. KOLLE HANDLE

Loop end.

Handle length: 43 mm. Loop Ø: 3 mm.

Part No.

**1000013** Made of Ni-Cr.

**1000014** Made of Platinum iridium (10% iridium).

#### 2. NEEDLE FOR KOLLE HANDLE

Made of Ni-Cr.

Part No. **1000802**

#### 3. LANCET TIP FOR KOLLE HANDLE

Made of Ni-Cr. Total long: 50 mm.

Part No. **1000803**




### STAND

Support rack for Kolle handles dissection needles made of hardened methacrylate.

Capacity: 6 places.

Part No. **1000903**



### STERILE CULTURE LOOPS

Disposable flexible polystyrene culture handle, one loop and one pointed end.

Part No.	Capacity	Colour	internal loop Ø	In packets of
<b>1000866</b>	1 µl	green	1 mm	20 units
<b>1000867</b>	10 µl	blue	4 mm	20 units



### INFRARED TIP STERILIZER "STERILBIO"

#### Applications

Sterilization of culture tips, needles and small instruments. Incinerates organic material at 900 °C for 5/8 seconds.

#### Features

Made from ceramic elements with a quality quartz tube interior.

Does not need purged oxygen and can be used in anaerobic environments.

Eliminates problems of aerosols and bacterial contamination.

Support at the side for 4 Kolle handles, needles, lancet or dissecting forceps etc.

#### Model

Part No.	Height/Width/Depth (exterior) cm			Ø/Length (tube interior) cm		Power W	Weight Kg
<b>3000788</b>	20	13	18	1.4	15	110	1.5