

**BOTTLES AND FLASKS: Plastic flasks**

**Plastic flasks**

Squared plastic flasks, made of high density polyethylene (HDPE), to store laboratory reagents, without cap, with GL45 thread which allows to thread dispensers and burettes. Due to its great sturdiness, prevents collapse with alcohols. Because of its base, they are very stable even if empty. Fit perfectly on automatic titration supports. Volume marks allow estimation of the remaining liquid in flask.

Cap. (ml)	Thread GL	Pack (u.)	Art. No.
100	45	20	FAB-0PI150
250	45	20	FAB-0PI250
500	54	20	FAB-0PI500
1000	54	20	FAB-PI1000
1.500	80	20	FAB-PI1500
2.500	80	20	FAB-PI2500


**Plastic bottles**

Plastic bottle made of high density polyethylene (HDPE) suitable for storing reagents, without cap.

Cap. (ml)	Thread GL	Pack (u.)	Art. No.
100	32	20	FAB-RI0100
250	32	20	FAB-RI0250
500	45	20	FAB-RI0500
1000	45	20	FAB-RI1000
2.500	45	20	FABRI2500N
4000	45	20	FAB-RI4000

**Economic package without cap**

Capacity (ml)	Thread PILFER	Pack (u.)	Art. No.
500	28	20	FAB-BAS500
1000	28	20	FAB-BA1000


**Caps for plastic flasks**

Caps for plastic flasks.

Thread GL	Pack (u.)	Art. No.
45	20	FAB-TA1337
54	20	FAB-TAPS54
80	20	FAB-TAPS80

**Caps for plastic bottles**

Thread GL	Pack (u.)	Art. No.
GL 45	20	FAB-TA1639
Pilfer 28 (for economic package)	20	FABTAP125B

Visit our new web [www.scharlab.com](http://www.scharlab.com)