

Interscience 1 2

BagMixer® 400. Range of 400ml lab blenders



World best-seller, the BagMixer® 400 is a powerful and easy-to-use lab blender. In use with a filter-bag, there is no risk of cross-contamination: the blending is accurate, reliable and quick. Adapted to any kind of application, it is the perfect tool for lab analysis.

Applications

- Food, pharmaceutical, clinical, veterinary and environmental analyses

BagMixer® 400 P The basic essentials

- Fixed speed (8 strokes/sec.)
- Variable blending time (30 - 210 sec. or open-running)
- Also available in CoolMix® version (refrigerated / heated door, from -15°C to 80°C)

BagMixer® 400 VW Comfort and precision

- Window door
- Variable speed (6 - 9 strokes/sec.)
- Variable blending time (10 - 360 sec. or open-running)
- Adjustable blending strength (adjustable paddles)
- Security drip tray

BagMixer® 400 W Performance

- Window door
- Fixed speed (8 strokes/sec.)
- Variable blending time (30 - 210 sec. or open-running)
- Security drip tray

BagMixer® 400 CC® Click and Clean®

- Window door
- Multifunction digital screen
- Variable speed (3 - 12 strokes/sec.)
- Variable blending time (1 sec - 1 h or open-running/ count-down)
- Adjustable blending strength with on-screen indication
- Click & Clean® System (removable paddles)
- Security drip tray

Order No.	Description
603-000004 1	BagMixer® 400P
603-000007 2	BagMixer® 400W
603-023230	BagMixer® 400VW

Technical specifications

BagMixer® 400	P	W	VW	CC®
Useful volume:	50 - 400 ml	50 - 400 ml	50 - 400 ml	50 - 400 ml
Power:	220 / 110 V - 50 / 60 Hz	220 / 110 V - 50 / 60 Hz	220 / 110 V - 50 / 60 Hz	220 / 110 V - 50 / 60 Hz
Size (L x l x h)	39 x 26 x 29 cm	39 x 26 x 29 cm	42 x 26 x 29 cm	42 x 26 x 26,5 cm
Weight:	16,5 kg	16,5 kg	16,5 kg	16,5 kg
Reference:	603-000004	603-000007	603-023230	603-024230

Interscience 3

BagMixer® Click & Clean®



The Click & Clean® system offers great user-comfort and unequalled blending quality. The blending chamber is 100% accessible for complete cleaning, due to the removable and autoclavable paddles.

Order No.	Description
603-024230	BagMixer® 400CC®. Window door

Click & Clean® system

The Click and Clean® system allows integral access to the blending chamber: in 1 click, remove the paddles, autoclave them, clean the blender.



Patented

Adjustable blending strength

Adjust the blending chamber according to the nature and size of your sample.



Patented

Easy control

- Blending speed
- Blending strength
- Blending time

