

PIPETTES: Pipetting filler

“FLIP” Safety pipette filler
DEUTSCH & NEUMANN

Safety pipette filler, natural rubber.
 Only two operating points
 Evacuate via the new automatic valve
 Ability to clean inside of bulb by removing patented valve

Colour	Pack (u.)	Art. No.
Red	1	0004710001
Black	1	0004710011


Rubber bulbs for smaller pipettes
DEUTSCH & NEUMANN

Description	Pack (u.)	Art. No.
Natural rubber	100	4600002DAN
Silicone rubber	100	4600003DAN


Safety pipette fillers
DEUTSCH & NEUMANN

Natural rubber pipette fillers.

Model	Colour	Pack (u.)	Art. No.
Standard, short neck	Red	1	0004700001
Universal, large neck	Red	1	0004700002
Universal, with adapter	Red	1	425-000200

For further information, please contact scharlab@scharlab.com