



Histology paraffin section mounting bath "Termofin"

ADJUSTABLE TEMPERATURE FROM 40 °C TO 80 °C.
STABILITY ± 1 °C.

APPLICATIONS

Suitable for the fusion and preservation of paraffin wax in a liquid state that enables the user to see the floated sections.

FEATURES

Aluminium inner tank and lid with a P.T.F.E. matt black coating. Thermometer location hole in the lid to read temperature.

CONTROL PANEL

Mains power ON/OFF switch. Hydraulic temperature thermostat regulator synchronised with an indicator lamp.

MODEL

Part No.	Capacity litres	\emptyset (usable) cm	Height cm	\emptyset (exterior) cm	Height cm	Power W	Weight Kg
3000459	4	25	10	34	14	430	2.3

Comes complete with lid and thermometer.

P.T.F.E. coated



ACCESSORY

Thermometer support (can only be used without the lid).

Part No. **6004591**

NOTE: Wax dispenser (see page 217).



Bank of oil and water baths "Precisbat"

ADJUSTABLE TEMPERATURES FROM AMBIENT +5 °C UP TO 200 °C, STABILITY ± 2 °C.
INDEPENDENT TEMPERATURE CONTROL FOR EACH POSITION.

FEATURES

Inner vessel made of AISI 304 stainless steel with the external case made of AISI 304 stainless steel.

The heater element is embedded in stainless steel INCOLOY alloy, that is corrosion and heat resistant.

Each position has a lid with concentric reduction rings of; 105, 80, 60 and 40 mm in \emptyset .

An adjustable support frame is mounted at the back of the unit to facilitate the mounting and support of flasks, condensers and other apparatus, etc.

Capacity: 2 litres per position.

CONTROL PANEL

Independent control for each position.

Mains ON/OFF with luminous "ON" indicator lamp for each position.

Hydraulic thermostat temperature control synchronised with the "ON" indicator lamp.



MODELS

Part No.	Temperature maximum °C	Number of places	Height / Width / Depth (usable each position) cm			Height / Width / Depth (exterior) cm			Power W	Weight Kg
6001462	200	2	12	14	12	19	35	29	1300	5.5
6001472	200	4	12	14	12	19	67	29	2560	10
6001482	200	6	12	14	12	19	100	29	3800	15

ACCESSORIES



Perforated shelf (One per tank) made of AISI 304 stainless steel. Part No. **6000011**



Clamp for supporting flasks and condensers. Part No. **7000403**



Thermometer up to 220 °C. Part No. **1000873**



Bosshead for supporting clamps. Part No. **1000069**