



### Incubator CULTURA M

Small incubator for Petri dishes, dipslides, etc...  
 Temperature range: 25 - 45°C.  
 Temperature stability: approx. 1°C.  
 Materials: PC & ABS. Door made of transparent PC.  
 Thermometer in special tube (up to 60°C).  
 Includes a removable rack for 18 tubes.  
 External dimensions (WxHxD): 310x155x168mm.  
 Weight: 1.1kg.

Description	Pack (u.)	Art. No.
CULTURA M incubator	1	000070700R



### Jar for anaerobic incubation

Jar of 2,5 liters, for 12 plates with basket. Safety seal.

Description	Pack (u.)	Art. No.
Jar for anaerobic incubation	1	064-AG0P25

### Anaerobiosis generator sachets

Do not require a catalyst.

Description	Pack (u.)	Art. No.
Anaerobic atmosphere generator sachets	10	230-096124
Microaerophilic atmosphere generator sachets	10	230-096125
Sachets for capnophilic bacteria	10	230-096126
Anaerobiosis Indicators	50	230-096118



### Mini incubator EM4 for slides

Small, lightweight unit, with a transparent cover and handle to ease extraction. Includes a rack of 23 dip-slides. Includes thermometer. Suitable for incubation of dip-slides, Petri dishes, RODAC plates, tubes, etc. Adjustable temperature from ambient +5°C up to +44°C. Dimensions (WxHxD): 23x18x23cm. Weight: 1,5Kg.

Description	Pack (u.)	Art. No.
Mini incubator EM4	1	114-000020

For further information, please contact [scharlab@scharlab.com](mailto:scharlab@scharlab.com)