

# Nagy méretű kerámia lapos fűtőlap 30 x 50cm analóg

## / BS-00CB5000



Ceramic hotplate surface

Large plate area ideal for heating multiple vessels

Accommodates up to 30 x 100ml beakers

This large rectangular shaped hotplate is ideal for heating many smaller vessels simultaneously, e.g. in educational use where lots of students require their samples heated at the same time.

An easy to use knob located on the front of the hotplate controls the top plate temperature. The robust side arms make the unit very easy to carry.

The large surface area may stay hot for a long time after use, so for maximum safety, a bright red hot warning light will continue to flash until the hotplate is cool. There are two models to choose from:

Model CB500 has a glass ceramic top plate. It is easy to clean in the event of spillage and, due to unique thermal properties, can be heated to very high temperatures giving very fast heat up times.

Full product specification.

	CB500
Plate material	Glass ceramic
Plate dimensions, mm	300 x 500
Heated area, mm	250 x 450

Heater power, W	2250
Max. plate temp. °C	375
Overall dimensions, (w x d x h)	520 x 360 x 130
Net weight, kg	12
Electricity supply	230V, 50-60Hz