

Bandelin

1 2

SONOREX DIGITEC Modern High-Power Ultrasonic Baths

- Cleaning of technical glassware like burettes, pipettes, petri dishes and laboratory flasks.
- Disinfection and cleaning at the same time.
- Degassing of beer samples for analysis of alcohol contents, original worth, colour, pH value.
- Degassing of food samples from cans for analysis of stannous contents.
- Extraction of quaternary ammonium compounds (QAC) of wood.
- Extraction of herbs samples for determination of aflatoxines (causing mold decay on food).
- Extraction of soil samples for determination of hydrocarbons.
- Test method for freeze-thaw resistance of concrete: CDF test through sonication, loosely adhering scaled particles are removed from surface.

1



2



More details!
 Use of insert gasket.
 Please note: Highly recommended to avoid damage of your bath.

NEW!

Type	Capacity	Inside dims. (L x W x D)		Outside dimensions (L x W x H)		Drain ball cock	HF-peak output W / Per.	HF peak W	Heating power W	Weight empty kg	PU	Order No.
		l	mm	mm	mm							
DT 31	0.9	190 x 85 x 60	205 x 100 x 170				240	30	-	1.8	1	0000003200
DT 31 H	0.9	190 x 85 x 60	205 x 100 x 170				240	30	70	1.9	1	0000003220
DT 52	1.8	150 x 140 x 100	175 x 165 x 230				240	60	-	2.6	1	0000003205
DT 52 H	1.8	150 x 140 x 100	175 x 165 x 230				240	60	140	2.9	1	0000003225
DT 100	3	240 x 140 x 100	260 x 160 x 250				320	80	-	3.4	1	0000003210
DT 100 H	3	240 x 140 x 100	260 x 160 x 250				320	80	140	3.6	1	0000003230
DT 102 H	3	240 x 140 x 100	260 x 160 x 250	G			480	120	140	4.3	1	0000003235
DT 103 H	4	240 x 140 x 150	260 x 160 x 310	G			560	140	200	4.6	1	0000003201
DT 106	5.6	∅ 240 x 130	∅ 265 x 270	G			480	120	-	5.5	1	0000003270
DT 156	6	500 x 140 x 100	530 x 165 x 245	G			640	160	-	6.1	1	0000003275
DT 156 BH	9	500 x 140 x 150	530 x 165 x 300	G			860	215	600	7.3	1	0000003221
DT 255	5.5	300 x 150 x 150	325 x 175 x 295	G			640	160	-	5.2	1	0000003215
DT 255 H	5.5	300 x 150 x 150	325 x 175 x 295	G			640	160	280	5.3	1	0000003240
DT 510	9.7	300 x 240 x 150	325 x 265 x 305	G			640	160	-	7.0	1	0000003245
DT 510 H	9.7	300 x 240 x 150	325 x 265 x 305	G			640	160	400	7.6	1	0000003206
DT 512 H	13	300 x 240 x 200	325 x 265 x 350	G			860	215	400	8.0	1	0000003226
DT 514	13.5	325 x 300 x 150	355 x 325 x 305	G			860	215	-	8.2	1	0000003250
DT 514 H	13.5	325 x 300 x 150	355 x 325 x 305	G			860	215	600	8.8	1	0000003211
DT 514 BH	18.7	325 x 300 x 200	355 x 325 x 385	G			860	215	600	9.8	1	0000003216
DT 1028	28	500 x 300 x 200	535 x 325 x 400	G			1200	300	-	14.3	1	0000003255
DT 1028 C	45	500 x 300 x 300	540 x 340 x 500	G			2000	500	-	24.6	1	0000003231
DT 1028 H	28	500 x 300 x 200	535 x 325 x 400	G			1200	300	1300	14.7	1	0000003231
DT 1028 CH	45	500 x 300 x 300	540 x 340 x 500	G			1200	300	1450	23.7	1	0000003266