



Heating cleaning bath “Ultrasons-HD”

DIGITAL ELECTRONIC REGULATION OF TEMPERATURE AND TIME
FOR TEMPERATURES FROM AMBIENT +5°C TILL 80°C

Optimum cleaning efficiency



It comes complete with a universal basket, lid with handle and drain device.

Bath Part. No. 3000865 (5 Litres capacity) is also supplied with basket for small items.

FEATURES

Double body container, completely made of stainless steel and with waste drain device.

CONTROL PANEL

Illuminated main switch.

Start-stop push button.

Increase/decrease time push button in minutes (from 0' to 60').

Increase/decrease temperature push button in °C.

MODEL

Part No.	Capacity litres	Height/Width/Depth (usable) cm			Height/Width/Depth (exterior) cm			Power Generator W	Power Heater W	Weight Kg
3000865	5	15	24	14	27	31	18	120	75	3,8
3000866	6	15	30	15	27	35	18	180	75	5,8
3000867	10	15	30	24	35	38	28	240	225	8,2
3000868	22	15	50	30	27	58	35	600	375	15
3000869	40	20	50	40	40	65	50	900	1000	25

It comes complete with a universal basket, small basket for small items, lid with handle and drain device.



Heated Ultrasonic baths “Ultrasons-H”

AJUSTABLE TEMPERATURE FROM AMBIENT +5 °C HASTA 75 °C.

APPLICATIONS

For long pieces cleaning (especially for pipettes) .



Note: For cleaning long objects whose length may exceed that of the tank, such as tubes, pipettes, etc, it has been proved that they are submitted to the same cavitation effect by introducing them in a test tube filled with liquid and placed vertically in the cleaner.

FEATURES

Evenly distributed heating element, adhered to the tank for maximum heat transfer.

CONTROL PANEL

Adjustable thermostat heater.

Heater “on” indicator lamp.

Timer switch from 0-15 minutes or continuous operation.

In operation on indicator lamp.

MODEL

Part No.	Capacity litres	Height/Width/Depth (usable) cm			Height/Width/Depth (exterior) cm			Power Generator W	Power Heater W	Weight Kg
3000839	9	15	50	14	30	56	20	200	800	11

ACCESSORIES

Basket, Height 12/Width 46/Depth 12 cm. Part No **6005144**

3 Lid with concentric reduction ring. Part No **6005134**

Lid with handle. Part No **6005124**

See page: 91



Ultrasonic pipette cleaner “Ultrasons-P”

FEATURES

Designed for automatically cleaning and rinsing pipettes, mounted inside a basket, pipettes are ultrasonically cleaned and rinsed.

Minimum rinsing time by using a water siphon action through the inside of the pipettes.

Metallic areas in contact with liquid are made of AISI 316 stainless steel.

The outer casing is made of AISI 304 stainless steel.

Cylindrical vessel in transparent acrylic, with lid to avoid splashing.

Connector(1”) water inlet and siphon-drain tube.

Basket with handle in polypropylene for the easy handling and carrying of the pipettes.

CONTROL PANEL

Timer from 0-15 minutes or continuous operation.

Indicator lamp showing when the bath is in operation.

MODEL

Part No.	Capacity litres	Ø / Height (usable) cm	Height / Width/ Depth (exterior) cm			Power W	Weight Kg
3000910	10	15 64	84	34	40	250	15