

# WIZARD

Vortex Mixer

## The Next Generation Vortex Mixer!

“NO PRESS, NO STRESS”

Stirring Line

A new solution that will optimise and simplify your daily tasks

### No Stress

Drastic reduction in repetitive strain

### Automatic Operation

No pressure required  
Proximity detector

### Exceptional Stability on the Bench

Low profile and special materials for excellent convenience and robustness

### Ergonomic and Highly Innovative Design

Awarded from the prestigious American publisher Strategic Directions International, Inc. (SDi), through its newsletter Instrument Business Outlook (IBO) for excellence in design that makes it unique

**Vortex “Wizard” is a Patented Solution**



*Industrial Design Award*

## Additional benefits of choosing the VELP WIZARD Vortex Mixer

**Optimized Performance**  
**Quiet Operation**  
**Usable on Many Surfaces**  
**Supreme Inherent Safety**  
**Versatile Adaptability**

0-3000 rpm + soft start + controlled ramping  
Low noise motor plus insulation  
3-point anti-skid feet  
IP42 – highest safety of any mixer on the market  
Wide range of accessories

Technical Data	Description
Construction material:	Zinc alloy and technopolymer
Support system:	3-point anti-skid feet
Type of movement:	Orbital
Orbital diameter:	4.5 mm
Operational mode:	Sensor / Continuous
Speed setting:	Analogue
Speed control:	Electronic
Stirring speed:	0 – 3000 rpm
Electrical protection degree (CEI EN 60529) :	IP 42
Power:	15 W
Power supply:	100 – 240 V / 50 – 60 Hz
Weight:	2.2 Kg
Dimension (W x H x D) :	180 x 70 x 220 mm

Ordering Information Code N°	Description
F202A0175	Vortex Mixer WIZARD
A00000012	Foam stand for n° 19 Eppendorf® type 1.5 ml microvials
A00000013	Customizable soft foam top
A00000014	Foam stand for n° 5 test tubes Ø 16 mm
A00000015	Foam stand for microtiter
A00000016	Small rubber supporting plate Ø 50 mm
A00000019	Foam stand for n° 4 test tubes Ø 29 mm
10003083	US plug
10003084	UK plug
10003085	Australian plug



A00000012



A00000013



A00000014



A00000019



A00000015




A00000016

Also available: **WX** model - Basic version with the “Sensor” mode only (code N° F202A0230).

For further information please contact [inse@velp.it](mailto:inse@velp.it).

We reserve the right to make technical alterations  
We do not assume liability for errors in printing, typing or transmission



VELP Scientifica srl  
via Stazione 16  
20040 Usmate (Milano) Italy  
Tel +39 039 628811  
Fax +39 039 6288120  
[inse@velp.it](mailto:inse@velp.it)  
[www.velp.com](http://www.velp.com)