



PP tube with screw cap

Test tubes with screw caps

KARTELL

Test tubes with quick closure screw caps (PE). Leak-proof, and suitable for centrifugation. The design of the cap avoids any contamination due to drops of liquid on the inner surface.

Description	Bottom	Cap. (ml)	Dim. ØxLength (mm)	Material	Tubes /box	Pack (u.)	Art. No.
PP tube with screw cap	Round	15	16x100	PP	1.000	1x1.000	425-088303
PP tube with screw cap. Sterile	Round	15	16x100	PP	100	1x100	425-088304
Tube with cap	Conical	50	40x120	PP	100	1x100	425-084002
Microtube cap for Cobas-Bio®	Conical	0,7	7,6x35,5	HDPE	1.000	1x1.000	4252940/00
Microtube cap for Cobas-Bio®	Conical	0,7	7,6x35,5	HDPE	1.000	1x1.000	4252940/04



PS tubes

Disposable tubes

KARTELL

Made of polystyrene or polypropylene round or conical tubes, no rim, for laboratory use. Can be used in centrifuge.

Description	Bottom	Cap. (ml)	Dim. ØxLength (mm)	Material	Tubes/box	Pack (u.)	Art. No.
PS test tubes. radioimmunoassay (RIA)	Round	3	11x70	PS	1.000	6x1.000	425-088310
PS test tubes	Round	5	12x75	PS	1.000	1x1.000	425-088307
PS test tubes	Conical	10	16x105	PS	2.000	1x2.000	425-088301
PS test tubes	Round	10	16x100	PS	2.000	1x2.000	425-088302
PS test tubes	Round	15	16x150	PS	1.000	1x1.000	425-088320
PS test tubes. Coagulometer	Round	11,5	11,5x55	PS	1.000	5x1.000	425-088308



Tube with screw cap and skirt

Test tubes with self standing

DELTALAB

Made of autoclavable (121°C) transparent polypropylene or polystyrene.

Dimensions: 15 x 102mm.

Green cap made of high density polyethylene. The external self standing allows the tubes to remain free-standing. Supplied screwed.

The sterile model (ethylene oxide) is supplied individually in flow-pack bag, with indications of batch, expiration date, etc...

Resistance to centrifugation: PS: 7.500g. PP: 15.000 g.

Attention: for autoclaving, the cap should be loose on the thread and not tightly fitted.

Description	Bottom	Capacity (ml)	Dim. ØxLength (mm)	Material	Tubes/box	Pack (u.)	Art. No.
Tube with screw cap and self standing	Round	12	15x102	PS	250	6x250	027-301403
Tube with screw cap and self standing	Round	12	15x102	PP	250	6x250	027-401403
Tube with screw cap and self standing. Sterile, unit bag	Round	12	15x102	PS	250	6x250	027-301402
Tube with screw cap and self standing. Sterile, unit bag	Round	12	15x102	PP	250	6x250	027-401402

For further information, please contact scharlab@scharlab.com