

## Specifikáció



### Termék kódja , megnevezése:

OCT0006600	Habzágátló tablettá (Antifoam) 0,03g Na <sub>2</sub> SO <sub>4</sub> - 0,97g S - szilikon 1000db/doboz
------------	---

### Reagents tablets

Catalysts and antifoam for accelerating, simplifying and standardizing oxidative mineralizations in nitrogen analyses according to Kjeldahl. Catalysts and antifoam are packaged in pre-dosed tablets; this helps prevent possible errors in the dosage of chemical substances and speeds up the analytical process considerably. They can easily be used in agricultural chemistry, clinical chemistry and petrochemistry, and for analyzing foods and animal feeds.