



PP tube with screw cap

### Test tubes with screw caps

#### KARTELL

Test tubes with quick closure screw caps (PE). Leak-proof, and suitable for centrifugation. The design of the cap avoids any contamination due to drops of liquid on the inner surface.

Description	Bottom	Cap. (ml)	Dim. ØxLength (mm)	Material	Tubes /box	Pack (u.)	Art. No.
PP tube with screw cap	Round	15	16x100	PP	1.000	1x1.000	425-088303
PP tube with screw cap. Sterile	Round	15	16x100	PP	100	1x100	425-088304
Tube with cap	Conical	50	40x120	PP	100	1x100	425-084002
Microtube cap for Cobas-Bio®	Conical	0,7	7,6x35,5	HDPE	1.000	1x1.000	4252940/00
Microtube cap for Cobas-Bio®	Conical	0,7	7,6x35,5	HDPE	1.000	1x1.000	4252940/04



PS tubes

### Disposable tubes

#### KARTELL

Made of polystyrene or polypropylene round or conical tubes, no rim, for laboratory use. Can be used in centrifuge.

Description	Bottom	Cap. (ml)	Dim. ØxLength (mm)	Material	Tubes/box	Pack (u.)	Art. No.
PS test tubes. radioimmunoassay (RIA)	Round	3	11x70	PS	1.000	6x1.000	425-088310
PS test tubes	Round	5	12x75	PS	1.000	1x1.000	425-088307
PS test tubes	Conical	10	16x105	PS	2.000	1x2.000	425-088301
PS test tubes	Round	10	16x100	PS	2.000	1x2.000	425-088302
PS test tubes	Round	15	16x150	PS	1.000	1x1.000	425-088320
PS test tubes. Coagulometer	Round	11,5	11,5x55	PS	1.000	5x1.000	425-088308



Tube with screw cap and skirt

### Test tubes with self standing

#### DELTALAB

Made of autoclavable (121°C) transparent polypropylene or polystyrene.

Dimensions: 15 x 102mm.

Green cap made of high density polyethylene. The external self standing allows the tubes to remain free-standing. Supplied screwed.

The sterile model (ethylene oxide) is supplied individually in flow-pack bag, with indications of batch, expiration date, etc...

Resistance to centrifugation: PS: 7.500g. PP: 15.000 g.

Attention: for autoclaving, the cap should be loose on the thread and not tightly fitted.

Description	Bottom	Capacity (ml)	Dim. ØxLength (mm)	Material	Tubes/box	Pack (u.)	Art. No.
Tube with screw cap and self standing	Round	12	15x102	PS	250	6x250	027-301403
Tube with screw cap and self standing	Round	12	15x102	PP	250	6x250	027-401403
Tube with screw cap and self standing. Sterile, unit bag	Round	12	15x102	PS	250	6x250	027-301402
Tube with screw cap and self standing. Sterile, unit bag	Round	12	15x102	PP	250	6x250	027-401402

For further information, please contact [scharlab@scharlab.com](mailto:scharlab@scharlab.com)