



Test tubes

Conical test tubes DELTALAB

Manufactured in polypropylene and polystyrene. Art. No.s 301200, 401200 and 401204 have internal graduations. Art. No.s 301213 and 301212 are flared at their top for easy pouring and have a special sedimentation base. Art. No.s 301200, 401200 and 301202 have an external rim. All models are suitable for centrifugation.
 PS tubes: 7.500g; PP tubes: 15.000g.
 Except Art. No. 301212 and 301213 to 3.000g.

Description	Bottom	Cap. (ml)	Dim. ØxLength (mm)	Material	Tubes/box	Pack (u.)	Art. No.
PS tube	Conical	12	16x102	PS	500	5x500	027-301200
PS tube	Conical	12	16x100	PS	500	5x500	027-301201
PS tube	Conical	12	16x102	PS	500	5x500	027-301202
PS tube	Conical	12	17x105	PS	1.500	1x1.500	027-301212
PS tube	Conical	12	16x105	PS	1.500	1x1.500	027-301213
PS tube	Conical	12	16x102	PP	500	2x500	027-401200
PS tube	Conical	12	16x100	PP	500	5x500	027-401201
PS tube	Conical	12	16x100	PP	500	5x500	027-401204

Round bottom polystyrene tubes DELTALAB

Tubes are made in one piece to ensure uniformity and dimensional accuracy from tube to tube. Art. No. 300800 is graduated at 1, 2,5 and 4,0ml; Art. No. 300400 is graduated at 1, 2,5 and 5,0ml. Suitable for centrifugation up to 7.500g. All tubes are suitable for laboratory use. Excepting Art. No. 300100, which does not withstand centrifugation.



Description	Bottom	Cap. (ml)	Dim. ØxLength (mm)	Material	Tubes/box	Pack (u.)	Art. No.
PS Tube. Without leveling	Round	5	12x75	PS	1.000	4x1.000	00300800.1
PS Tube. Thick wall	Round	5	12x75	PS	1.000	4x1.000	00300800.2
PS Tube	Round	4	11x70	PS	1.000	3x1000	027-300300
PS Tube	Round	6	12x88	PS	1.000	4x1.000	027-300400
PS Tube	Round	3	11x55	PS	2.000	2000	027-300500
PS Tube	Round	5	12x75	PS	1.000	4x1.000	027-300800
PS Tube	Round	17	16x150	PS	1.000	1x1.000	027-300100
PS Tube	Round	5	13x75	PS	1.000	4x1.000	027-300700
PS Tube	Round	11	16x100	PS	500	5x500	027-300705
PS Tube	Round	10	16x95	PS	500	5x500	027-300900
PS Tube	Round	7	13x100	PS	500	6x500	027-301700
Capped tube	Round	6	13x75	PS	2.000	2.000	027-300702
Capped tube	Round	9	15x100	PS	500	500	027-300907