


**Conical tubes 15ml polypropylene
DELTALAB**

Description	Dim. Ø ext.xH (mm)	Resist. centrif. (g)	Sterile	Pack (u.)	Art. No.
15ml. Molded graduation and writing area	17x120	5.000	No	500	027-429910
15ml. Molded graduation and writing area	17x120	5.000	Sterile	5x100	027-429920
15ml. High resistance. Printed graduation. Free of DNase, RNase, endotoxins and metals	17x118	17.000	Sterile	10x50	027-409920

**Conical tubes 50ml polypropylene
GOSSELIN SAS**


Description	Dim. Ø ext.xH (mm)	Resist. centrif. (g)	Sterile	Pack (u.)	Art. No.
50ml. Conical base. Printed graduation and writing area	29,2x113,5	9.600	No	20x25	630TC50-01
50ml. Conical base. Printed graduation and writing area	29,2x113,5	9.600	Sterile	20x25	630TC50-04
50ml. Self standing. Printed graduation and writing area	29,2x113,5	9.600	No	20x25	00TU50-012
50ml. Self standing. Printed graduation and writing area	29,2x113,5	9.600	Sterile	20x25	00TU50-013


**Conical tubes 50ml polypropylene
DELTALAB**

Description	Dim. Ø ext.xH (mm)	Resist. centrif. (g)	Sterile	Pack (u.)	Art. No.
50ml. Conical base. Molded graduation	29,5x118,5	12.000	No	500	027-429900
50ml. Conical base. Molded graduation. Unit bag	29,5x118,5	7.000	Sterile	500	027-429926
50ml. Self standing. Printed graduation. High resistance. Free of DNase, RNase, endotoxins and metals	29,6x114,6	20.000	Sterile	10x50	027-409922
50ml. Self standing. Molded graduation	29,5x118,5	12.000	No	500	027-429901
50ml. Self standing. Molded graduation. Unit bag	29,5x118,5	7.000	Sterile	500	027-429927
50ml. High resistance. Printed graduation. Free of DNase, RNase, endotoxins and metals	29,6x114,6	10.000	Sterile	10x50	027-409923

 Visit our new web www.scharlab.com