



Scharlab s.l.
Tel. int.: +34-93-7151811
Email: lab@scharlau.com

CERTIFICATE OF ANALYSIS

Product: **Dichloromethane**, spectroscopy grade, stabilized with ethanol, Spectrosol® Batch 10991202
CL0333 Test date 10/09/2010
Shelf life 9/2013

Analysis	Batch value	Guaranteed value
assay (G.C.)	99.8 %	min. 99,8 %
identity (IR-spectrum)	passes test	passes test
density (20°C)	1.325	1,323 - 1,325
colour (Hazen)	8	max. 10
acidity	0.0001 meq/g	max. 0,0005 meq/g
alkalinity	< 0,0002 meq/g	max. 0,0002 meq/g
ethanol (G.C.)	0.2 %	max. 0,2 %
non-volatile matter	0.0001 %	max. 0,0002 %
water (K.F.)	0.005 %	max. 0,01 %
minimum transmission /max. absorbance in a 1,0 cm cell at		
wavelength:	T (%) A (AU)	T (%) A (AU)
230 nm	10 % 1,000	10 % 1,000
240 nm	70 % 0,155	70 % 0,155
250 nm	95 % 0,022	95 % 0,022
255 nm	98 % 0,009	98 % 0,009

This certificate does not release the user from the reception control .

This certificate is an electronic copy of the certificate available in our laboratory, and does not require signature

If you need further details, please call at our factory or contact our local distributor.

You can get a copy of any of our COA from our web site: www.scharlau.com