



**Scharlab S.L.**

Tel. int.: +34-93-7151811

Email: lab@scharlau.com

## CERTIFICATE OF ANALYSIS

Product: **Potassium hexacyanoferrate(II) trihydrate**, reagent Batch 12730902  
**P00248** grade, ACS, ISO, Reag. Ph Eur Test date 27/06/2011  
Shelf life 6/2016

Analysis	Batch value	Guaranteed value
assay (permanganometric)	100,51 %	99 - 102 %
insoluble in water	< 0,005 %	max. 0,005 %
chlorides (Cl)	< 0,01 %	max. 0,01 %
sulfates (SO <sub>4</sub> )	< 0,005 %	max. 0,005 %
cadmium (Cd)	< 0,0005 %	max. 0,0005 %
copper (Cu)	< 0,002 %	max. 0,002 %
lead (Pb)	< 0,002 %	max. 0,002 %
sodium (Na)	< 0,01 %	max. 0,01 %

This certificate does not release the user from the reception control .

This certificate is an electronic copy of the certificate available in our laboratory, and does not require signature

If you need further details, please call at our factory or contact our local distributor.

You can get a copy of any of our COA from our web site: [www.scharlau.com](http://www.scharlau.com)