



Scharlab s.l.
Tel. int.: +34-93-7151811
Email: lab@scharlau.com

CERTIFICATE OF ANALYSIS

Product: **Potassium sodium tartrate tetrahydrate**, reagent
P00355 grade, ACS, ISO

Batch 11794702
Test date 19/05/2010
Shelf life 5/2015

Analysis	Batch value	Guaranteed value
assay (alkalimetric)	100,4 %	99 - 102 %
insoluble in water	< 0,005 %	max. 0,005 %
pH (5 %, H ₂ O)	7,0	7,0 - 8,5
chlorides (Cl)	< 0,0005 %	max. 0,0005 %
phosphates (PO ₄)	< 0,001 %	max. 0,001 %
sulfates (SO ₄)	< 0,005 %	max. 0,005 %
total N	< 0,001 %	max. 0,001 %
ammonium (NH ₄)	< 0,002 %	max. 0,002 %
calcium (Ca)	< 0,004 %	max. 0,004 %
copper (Cu)	< 0,0005 %	max. 0,0005 %
heavy metals (as Pb)	< 0,0005 %	max. 0,0005 %
iron (Fe)	< 0,0005 %	max. 0,0005 %
lead (Pb)	< 0,0005 %	max. 0,0005 %
nickel (Ni)	< 0,0005 %	max. 0,0005 %
zinc (Zn)	< 0,0005 %	max. 0,0005 %

This certificate does not release the user from the reception control .

This certificate is an electronic copy of the certificate available in our laboratory, and does not require signature

If you need further details, please call at our factory or contact our local distributor.

You can get a copy of any of our COA from our web site: www.scharlau.com