

**Bandelin**

**SONOREX SUPER Compact Units**

1 2 3 4 5

PZT-transducers produce the high-powered ultrasonic intensity throughout the entire oscillating tank. Concentrated ultrasonic energy without irritating surface waves. Constant power and automatic frequency control always ensure optimum distribution of ultrasonic energy and reprodu-

cible results, regardless of the sizes of the pieces to be cleaned. Especially intense due to the advantageous frequency of 35 kHz.

**Features:**

- PZT-oscillating systems.
- HF-frequency 35 kHz.
- Timer 1 to 15 min and continuous operation.
- H-types with built-in heating, 30 to 80 °C thermostatically adjustable.
- Oscillating tank of stainless steel.

- Housing of stainless steel: drip proof.
- Drain cock included for easy emptying (from RK 102 H upward).
- Mains supply: 230 V, on request 115 V, 50 / 60 Hz.
- SONOREX SUPER units are RFI-proof and CE marked acc. to MPG.

**More details!**

Use of insert gasket.  
Please note: Highly recommended to avoid damage of your bath.



Type	Inside dimens. (L x W x D)	Capacity	Drain ball cock	HF-peak output W / Per.	HF peak W	Heating power W	Outside dimensions (L x W x H)	Weight	PU	Order No.	
	mm	l					mm	kg			
	RK 31	190 x 85 x 60	0.9	-	240	30	-	205 x 100 x 170	1.8	1	0000000329
<b>1</b>	RK 52	150 x 140 x 100	1.8	-	240	60	-	175 x 165 x 225	2.6	1	0000000311
	RK 52 H	150 x 140 x 100	1.8	-	240	60	140	175 x 165 x 225	2.9	1	0000000164
	RK 100	240 x 140 x 100	3.0	-	320	80	-	260 x 160 x 250	3.4	1	0000000301
<b>3</b>	RK 100 H	240 x 140 x 100	3.0	-	320	80	140	260 x 160 x 250	3.6	1	0000000312
	RK 100 SH	240 x 140 x 100	3.0	G 1 / 4	320	80	140	260 x 160 x 250	4.0	1	0000000192
<b>4</b>	RK 102 H	240 x 140 x 100	3.0	G 1 / 4	480	140	140	260 x 160 x 250	4.3	1	0000000303
	RK 103 H	240 x 140 x 150	4.0	G 1 / 4	560	140	200	260 x 160 x 310	4.6	1	0000000326
	RK 106	ø 240 x 130	5.6	G 1 / 4	480	120	-	ø 265 x 270	5.5	1	0000000306
	RK 156	500 x 140 x 100	6.0	G 1 / 4	640	160	-	530 x 165 x 245	6.1	1	0000000305
<b>2</b>	RK 156 BH	500 x 140 x 150	9.0	G 1 / 4	860	215	600	530 x 165 x 300	7.3	1	0000000646
	RK 255 H	300 x 150 x 150	5.5	G 1 / 4	640	160	280	325 x 165 x 305	5.3	1	0000000316
	RK 510	300 x 240 x 150	9.7	G 1 / 2	640	160	-	325 x 265 x 305	7.0	1	0000000327
	RK 510 H	300 x 240 x 150	9.7	G 1 / 2	640	160	400	325 x 265 x 305	7.6	1	0000000321
	RK 512 H	300 x 240 x 200	13.0	G 1 / 2	860	215	400	325 x 265 x 350	8.0	1	0000000795
	RK 514	325 x 300 x 150	13.5	G 1 / 2	860	215	-	355 x 325 x 305	8.2	1	0000000277
	RK 514 H	325 x 300 x 150	13.5	G 1 / 2	860	215	600	355 x 325 x 305	8.8	1	0000000207
	RK 514 BH	325 x 300 x 200	18.7	G 1 / 2	860	215	600	355 x 325 x 385	9.8	1	0000000263
	RK 1028	500 x 300 x 200	28.0	G 1 / 2	1200	300	-	535 x 325 x 400	14.3	1	0000000322
	RK 1028 H	500 x 300 x 200	28.0	G 1 / 2	1200	300	1300	535 x 325 x 400	14.7	1	0000000324
	RK 1028 C	500 x 300 x 300	45.0	G 1 / 2	2000	500	-	540 x 340 x 500	24.6	1	0000000661
	RK 1040	ø 500 x 195	38.0	G 1 / 2	1200	300	-	ø 540 x 500	20.5	1	0000000319