

## Narrow neck bottles series 301 LDPE - round, natural colour, with closure



Flexible, squeezable, with good shape-retaining properties, leakproof.

Also available with dropper cap or wash bottle closure.

### Application:

- \* dispensing \* adhesives \* sending samples
- \* transport \* storage

The packaging concept of Kautex for narrow, medium and wide neck laboratory bottles consists of consumer-optimized order units. Closures and bottles are packed together as one unit and can be ordered in an uncomplicated manner with one article number. This results in the following advantages for both dealers and end-users.


- \* one article number for bottle and closure
- \* better and clearer handling of the ordered goods
- \* dust-free and undamaged packaging
- \* easier and quicker to commission bottles and closures
- \* the end-user receives original packed merchandise

Nominal capacity ml	A mm	B mm	C mm	D mm	Mould no.	No. per case* pcs.	No. per packg. pcs.	Order no.
10	26	9,2	14	46,5	301-10	3200	100	301-770500
20	31,5	9,2	14	56,5	301-20	1800	100	301-770501
30	33,5	9,2	14	68,5	301-30	1400	100	301-770502
50	37,5	12,9	18	86	301-50	850	50	301-770503
100	47	12,9	18	105	301-100	480	20	301-770504
200	57,5	12,9	18	133	301-200	300	10	301-770505
250	61,3	19	25	140	301-250	210	10	301-770506
500	75	19	25	177	301-500	120	10	301-770507
1000	97,5	21,5	28	223,5	301-1000	56	4	301-770508
2000	120	21,5	28	264	301-2000	30	-	301-770509
3000	138	23,8	32	300	301-3000	20	-	301-770510
5000	162	30,8	40	355	301-5000	10	-	301-70511

\* Packaging unit / minimum ordering quantity

### Packaging concept:

- up to 100 ml: Closures are packed in a polybag which is included in the bottles polybag.
- ab 200 ml: Bottles are stretched in foil, corresponding number of closures in a separate polybag in the box.
- ab 2000 ml: Bottles packed loosely in the box with screwed-on closures.



①	<b>Dropper cap</b> with nozzle, natural colour
②	<b>Captive cap PVC</b> , red
③	<b>Wash bottle closure LDPE</b> , natural colour
④	<b>Screw closure LDPE</b> , natural colour (substitute, without picture)

### Closures

for bottle Nom. cap. ml	C mm	① Order no.	② Order no.	③ Order no.	④ Order no.
10/20/30	14	455-72102	419-73350	-	401-70603
50	18	455-72101	419-73350	-	401-70604
100	18	455-72101	419-73350	409-75991	401-70604
200	18	455-72101	419-73350	409-75992	401-70604
250	25	455-72100	419-73350	409-70618	401-70605
500	25	455-72100	419-73350	409-75993	401-70605
1000	28	455-72103	406-72103	409-70619	401-70606
2000	28	455-72103	406-72103	409-75994	401-70606
3000	32	-	-	-	401-70607

A = Diameter body B = Neck i.d. C = Cap size o.d. D = Height without closure